

L8 kit

TEST KIT TO PREVENT LEGIONELLA BACTERIA GROWTH

DESCRIPTION

Since the first reported cases of Legionnaires' Disease in 1976, outbreaks of this potentially fatal disease remains a threat. This has highlighted the need for frequent monitoring and testing of potential sources of infection.

The L8 KIT helps monitor the conditions that may allow legionella to grow. It has been assembled in-line with the new Health and Safety Commission (HSC) guidelines on legionella which require that specific, weekly tests be carried out on potential sources of infection. All the components needed for these tests are provided in the kit.

What is Legionnaires' Disease and where is it found?

Legionnaires' Disease is a type of pneumonia caused by a bacterium Legionella Pneumophila which is commonly found in natural and artificial water systems and can proliferate in poorly maintained wet cooling systems in a wide variety of workplaces (e.g. hospitals, office buildings, industrial premises and in the hotel and leisure industry).

Certain materials have been known to harbour the organism such as shower heads and rubber gaskets used in hot water systems. Infection may occur when fine water droplets containing bacterium are breathed in. It can also be found in cooling towers and evaporative condensers such as air conditioning plants and industrial cooling systems where typical temperatures may rise above 20°C.

Measures to control Legionnaires' Disease include proper disinfection and chlorination, thorough draining and regular inspection of all cooling towers and evaporative condensers.

Cisterns and pipework should be designed so that water is not allowed to stagnate or stand undisturbed for long periods.

An important part of reducing risks is a comprehensive monitoring programme. Weekly checks required by the HSC guidelines for cooling towers are:

- Oxidising biocide
- Conductivity
- pH
- Micro bacteriological activity

PRODUCT ORDERING INFORMATION

code SP 934L8

The L8 KIT includes:

TDS Sensor and calibration solution, comparator with disks and tablets for Bromine and Universal pH, test tube, comparator tube, brush, battery and 8 dipslides (2 supplied in box separately).



BLUTEC L8 KIT